

Supporting Information for “Hands-on Inquiry-Based Qualitative Identification of Metals in Coins Utilizing Atmospheric Pressure Chemical Ionization Mass Spectrometry”

Jonathan G. Moloney, Craig D. Campbell, Andrew F. Worrall\* and Malcolm I. Stewart




Chemistry Teaching Laboratory, South Parks Road, Oxford, U.K., OX1 3PS

\*Email: andrew.Worrall@chem.ox.ac.uk


## ***A117: Metal Chelate Complexes – Synthesis and Analysis***

## Introduction

- Acetylacetonate (acac) complexes have varied uses:  
e.g. catalysts for polymerisation;  $\text{Al}(\text{acac})_3$  is used as a molluscocide
- Traditionally prepared by replacing one ligand with acetylacetonate:







$\text{M}(\text{acac})_3$



$\text{Cu}(\text{acac})_2$

## Introduction

- Alternatively, it is possible to prepare the complexes electrochemically:



*Undivided  
cell*



## *Objectives*

- To investigate two different methods for the synthesis of a series of acetylacetonate complexes,  $M(\text{acac})_n$ :  
Salt metathesis and electrochemical synthesis
- To analyse the synthesised  $M(\text{acac})_n$  complexes by mass spectrometry
- To develop and investigate a method to identify metals present in coins

## ***Day 1 – Conventional Synthesis (Salt Metathesis)***

Each pair should prepare complexes in two categories:

- $\text{Fe}(\text{acac})_3$
- $\text{Cu}(\text{acac})_2 \cdot 2\text{H}_2\text{O}$
- $\text{Zn}(\text{acac})_2 \cdot \text{H}_2\text{O}$

These can be prepared using  $\text{MCl}_n$  or  $\text{M}(\text{OAc})_n$

$\text{NaOAc}$  or  $\text{NaOH}$  is used as the base to deprotonate acetylacetone  $\rightarrow \text{acac}^-$

Also prepare one of the following:

- $\text{V}(\text{acac})_3$  **[requires exclusion of air using a balloon of  $\text{N}_2$ ]**
- $\text{Cr}(\text{acac})_3$
- $\text{Mn}(\text{acac})_3$

Careful control of pH and/or adjustment of oxidation states required from most stable available metal salt

## *Considerations for Salt Metathesis*

Transition metal salts often have variable oxidation states

- Can be sensitive to oxidation in air/solution
- Can be sensitive to pH in aqueous solution, e.g. formation of insoluble  $M(OH)_n$

**V(acac)<sub>3</sub>:**  $V^{4+}$  oxidation state more stable – reduced to  $V^{3+}$  using Sn

**Cr(acac)<sub>3</sub>:**  $Cr^{3+}$  is sensitive to basic aqueous conditions

Urea is used to maintain weakly basic conditions:  $CO(NH_2)_2$  hydrolysed slowly to  $NH_3$  and  $CO_2$



$NH_3$  is used as a weak base to deprotonate acetylacetone

**Mn(acac)<sub>3</sub>:**  $Mn^{3+}$  salts often not available (aerobic stability)

Prepared by oxidation of  $Mn(acac)_2$  by  $MnO_4^-$

## *Electrochemical synthesis*

- You will choose one metal from Fe, Cu, and Zn to prepare  $M(\text{acac})_n$  electrochemically:



- NB: Although two electrodes are both the same metal, **oxidation** occurs only at the **anode (+)**, and **H<sub>2</sub> evolution** occurs at the **cathode (-)**
- Acetylacetone is mildly acidic ( $pK_a \sim 9$  in  $\text{H}_2\text{O}$ ); however, not sufficiently conductive
- Need an inert electrolyte to make the system significantly more conductive:  $[n\text{Bu}_4\text{N}] [\text{PF}_6]$

## *Safety*

### *Chemical safety*

- $\text{MX}_n$  compounds are corrosive/toxic
- Acetylacetone is flammable; methanol is flammable/toxic
- Keep electrolysis reactions in a fume-hood

### *Use of high voltage power supply*

- Do **not** allow the two electrodes to touch when power supply is on (sparks and ignites highly flammable organic solvent)
- Do **not** touch electrodes when power supply is turned on

### *Cleaning metal surfaces:*

- $\text{HCl(aq)}$  and  $\text{NaOH}$  are mildly corrosive





## *Summary*

### *Day 1:*

1. Prepare two  $M(\text{acac})_n$  complexes using conventional synthesis
2. Prepare one of  $\text{Fe}(\text{acac})_3$ ,  $\text{Zn}(\text{acac})_2$ , or  $\text{Cu}(\text{acac})_2$  electrochemically
3. Analyse  $M(\text{acac})_n$  by mass spectrometry

### *Day 2:*

Plan and investigate how to identify metals present in coins

**Complete the pro forma report for this practical and submit electronically to Canvas.**

***A117:  
Metal Chelate Complexes –  
Synthesis and Analysis  
Day 2***

## ***Mass Spectrometric Analysis***

- You will acquire the mass spectra of the synthesised  $M(\text{acac})_n$  complexes
- Mass spectrometry is a sensitive technique – our instrument can identify individual isotopes of each element to 0.1 Da
- APCI method: Soft ionisation technique. Generally detect solvated adducts, with possibly some fragmentation. Detect charged species of the forms:
- $[M+H]^+$ ,  $[M+H+H_2O]^+$ ,  $[M+H+CH_3OH]^+$ ,  $[M+H+CH_3COCH_3]^+$   
(and higher solvated adducts)



*i.e. [complex + H]<sup>+</sup>*




*i.e. [complex – acac]<sup>+</sup>*

Also possible: *Charged species* +  $H_2O$ ; + MeOH; + acetone

## ***Mass Spectrometric Analysis***

- You are provided with a spreadsheet of atomic masses and abundances of all naturally-occurring isotopes of the elements, “M”
- Use these data to determine the expected mass for the likely candidates that could be identified by mass spectrometry:




---

### *Sample preparation:*




Prepare a very dilute solution of your  $M(acac)_n$  in a vial in either MeOH or acetone – this will influence whether you may see +MeOH or +acetone adducts in the mass spectrum

## Mass Spectrometric Analysis

- For the synthesised  $M(\text{acac})_n$  you believe to have identified using the spreadsheet:
- Draw the structure for its charged species in ChemDraw
- Obtain its expected isotope composition  
(NB.  $^{13}\text{C}$  accounts for ~1.1% abundance of all carbon, which is not accounted for by your spreadsheet)

- View → Show Analysis Window

- e.g.  $[\text{Cu}(\text{acac})_2 + \text{H}]^+$  :



Major  $m/z$  species:

**262.2**

**264.0**

- Confirm the isotope pattern obtained is in agreement with the ChemDraw prediction