



FMRIB 3T MRI Scanning Procedure OxBAC-BIRAX Healthy Ageing (FMRIB 3T)

Study Information

Study Title: Healthy controls for the Oxford Health Centre: harmonising with UK Biobank

Study Group: Brain Health Centre / Neurology and Psychiatry

PI: [REDACTED]

Calpendo ID: 2021_015

Ethics Number: [REDACTED]

Introduction: Acquisition protocol for the studies: Oxford UK Biobank Aligned Controls (OxBAC) and British Israel Research Academic Exchange (BIRAX) study on Microstructural MRI in Healthy Ageing. The protocol includes sequences matched to the UK Biobank imaging study as well as specific sequences to study the brain microstructure, including structural connectome and cortical layers imaging.

General Info: Healthy participants aged 21 and over

Scanner: FMRIB 3T Prisma

Category: Category 1a - Healthy Subjects

Staffing levels: 1x scanner operator and 1x researcher

Session Time: 90 minutes

Incidental Findings Procedure: WIN SOP Dealing with Research Neuroimaging Incidental Findings (OHBA_014_V1, FMRIB_002_V5, Neuro_002_V5)

Researchers and their contact details

[REDACTED]

Equipment Required

Standard

32 Channel Head Coil
Regular mirror
Ear plugs
Immobilization sponges
Buzzer

Lighting and stimulus

Bore lights on
Room lights on
BOLD LCD screen on (pictures and fixation cross)


Ancillary

None


Protocol Sequences

FMRIB USER >> 2021_015 >> BIRAX >> Protocol v3

- 1. localiser_3plane_32ch** **00:14**
 - Auto-runs, no adjustment required
- 2. ep2d_diff_3scan_trace_p2_130** **00:43**
 - Axial orientation, angle to ACPC, straighten on coronal if required
 - Cover whole brain




- 3. TE_10_TE_20_TR30_...._PF78** **04:46**
 - Axial orientation, straighten on coronal if required
 - Copy COSG&SR from **ep2d_diff....**
 - Cover whole brain (move slice slab if necessary)



4. T1_p2_1mm_fov256_sag_TI_880 (biobank)

04:54


- Sagittal orientation, straighten on coronal & axial if required
- Cover whole brain, position to reduce phase wrap



5. t2_flair_sag_p2_1mm_FS_ellip_pf78

05:52


- Auto-copies COSG&SR from T1_p2_1mm.... (same FoV as T1)



6. TOF_3D_neck

00:42

- True axial, don't angle or straighten
- Position to cover from C2/C3 to lower cerebellum
- Once completed load into the viewer and find the best slice between vertebral loops




7. -- Manual shim x3 --

8. diff_PA_MPopt_MB3_3b0_lowflip

00:36

- Axial orientation, angle to ACPC, straighten on coronal if required
- Cover whole brain
- Open adjustments, click the 3D Shim tab and do 3x manual shims (ie Measure, Apply (top), Calculate, Apply (bottom))
- NB do not need to swap phase direction (already set in seq – special)
- Warn patient this scan, and the next, are very noisy and they may feel the table vibrate



9. diff_AP_MPopt_MB3_50b1000_50b2000_8b0_lowflip

06:32


- Auto-copies COSG&SR from diff_PA...
- NB do not need to swap phase direction (set in seq – special)

10. -- Copy diff Adj Vol --

11. to_ep2d_PCASL

07:17

- Axial orientation, angle to ACPC and adjust Sat region 1 to match
- Copy adj vol from diff_PA...
- Go to the Geometry Tab, select the Saturation Sub-Tab, choose Sat. region 1 and centre and angle to cover the scan volume
- Select Sat. region 2 and click the ... button, input the height chosen above, eg F76.8 and click enter in the centre of XX box
- Ensure the orientation is set to “transverse” If not, click on the three dots next to the orientation and reset it to “transverse”
- NB do not angle Sat region 2
- Resting state instructions: “Look at the cross for the next scan, blink normally and try not to fall asleep.”




12. -- Copy diff Adj Vol --

13. MB8_FMRI_fov210_2.4mm_resting

06:10

- Copy COSG&SR from #9 diff_PA_MPopt_MB3_3b)_lowflip...
- Copy Adj Vol from diff_PA...
- Cover whole brain
- Resting state instructions: “Look at the cross for the next scan, blink normally and try not to fall asleep.”




14. -- BIRAX --

15. dMRI_MB4_185dirs_d15D45_AP

11:10

- Axial orientation, angle to ACPC, straighten on coronal if required
- Cover cortex at expense of cerebellum



16. dMRI_MB4_6dirs_d15D45_APrev

00:44

- Auto copies COSG&SR from dMRI_MB4_185dirs_d15D45_AP
- Do not need to swap AP/PA


17. -- BIRAX (Optional) --

- The following scans take 12 minutes to run and can be skipped if running over time

18. IR-EPI_3iso_T150

00:12

- Axial orientation, angle to ACPC straighten on coronal as required
- Cover cortex at expense of cerebellum



19. IR-EPI_3iso_T1125

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

20. IR-EPI_3iso_T1225

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

21. IR-EPI_3iso_T1275

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

22. IR-EPI_3iso_T1325

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

23. IR-EPI_3iso_T1375

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

24. IR-EPI_3iso_T1410

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

25. IR-EPI_3iso_T1440

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

26. IR-EPI_3iso_T1460

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

27. IR-EPI_3iso_T1490

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

28. IR-EPI_3iso_T1510

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

29. IR-EPI_3iso_T1540

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

30. IR-EPI_3iso_T1570

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

31. IR-EPI_3iso_T1600

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

32. IR-EPI_3iso_T1620

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

33. IR-EPI_3iso_T1640

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150

34. IR-EPI_3iso_T1660

00:12

- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
35. IR-EPI_3iso_TI680	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
36. IR-EPI_3iso_TI700	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
37. IR-EPI_3iso_TI720	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
38. IR-EPI_3iso_TI740	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
39. IR-EPI_3iso_TI760	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
40. IR-EPI_3iso_TI780	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
41. IR-EPI_3iso_TI800	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
42. IR-EPI_3iso_TI820	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
43. IR-EPI_3iso_TI840	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
44. IR-EPI_3iso_TI860	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
45. IR-EPI_3iso_TI880	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
46. IR-EPI_3iso_TI910	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
47. IR-EPI_3iso_TI950	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
48. IR-EPI_3iso_TI990	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
49. IR-EPI_3iso_TI1010	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
50. IR-EPI_3iso_TI1050	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
51. IR-EPI_3iso_TI1100	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150]	
52. IR-EPI_3iso_TI1150	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
53. IR-EPI_3iso_TI1200	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
54. IR-EPI_3iso_TI1250	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
55. IR-EPI_3iso_TI1400	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
56. IR-EPI_3iso_TI1550	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
57. IR-EPI_3iso_TI1700	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	
58. IR-EPI_3iso_TI2000	00:12
- Auto copies COSG&SR from #16 IR-EPI_3iso_T150	

- | | |
|---|--------------|
| 59. IR-EPI_3iso_TI2100 | 00:12 |
| - Auto copies COSG&SR from #16 IR-EPI_3iso_T150 | |
| 60. IR-EPI_3iso_TI2500 | 00:12 |
| - Auto copies COSG&SR from #16 IR-EPI_3iso_T150 | |
| 61. IR-EPI_3iso_TI3000 | 00:12 |
| - Auto copies COSG&SR from #16 IR-EPI_3iso_T150 | |

62. -- Do retro-recon of SWI --

Post Scan

- Do retro-recon of SWI scan (see Appendix)
 - Researcher must confirm that the SWI retro-recon has been done before leaving
 - Researcher to return stimulus equipment to default state
-

Researcher Acknowledgement


Use of limited release sequences


The use of limited release sequences entails a risk that they will crash or otherwise slow down a session.

It is important to note that in the case where there is a delay caused by these sequences you may be required to cut other scans in order to finish within your allocated time.

Appendix – Retrospective Recon of SWI Data


1. The patient browser only shows the online recon scans




2. Shift+click the  button at the bottom of the screen to bring up the Retro Recon Window




3. The Retrospective Reconstruction window should now appear




4. Select the scan you want to re-reconstruct, then press “Expert edit”




5. The XBuilder window should now appear




6. Find the “ucCoilCombineMode” entry under the “MEAS” node, and change its value to “1”




7. Find the “ucUncomblImages” entry under the “MEAS” node, and change its value to “1”




8. Find the “ucSWI” entry under the “MEAS” node, and change its value to “0”



9. Save changes in XBuilder and close the window



10. Back in the Retro Recon window, press the "Start Recon" button to start reconstruction



11. In a few minutes the new images should appear in the Patient Browser

